

ELEVATORS

LUXATING ELEVATORS

Luxating elevators have thin, sharp blades that are used to cut periodontal ligaments. The luxating elevator is used by rocking and cutting vs. an elevator that lifts and pries. Luxating elevators reduce putting pressure on adjacent teeth.



Sharpening
Stone
Fine Grit
Ceramic

SSLUX

Straight
3mm

EL3S

Straight
5mm

EL5S

Curved
3mm

EL3C

Curved
5mm

EL5C

D	Wurzelheber
F	Elévateurs
I	Leve
E	Elevadores
P	Elevadores/Alavancas
J	エレベーター

DID YOU KNOW?

Luxating elevators are also available in a kit containing all of the above luxating elevators and a sharpening stone. The part code is ELKIT.